PHONAK | Government Services

Charger compatibility guide.

	Charger Go RIC E	Charger Go Virto I	Charger Go RICSphere I	Charger Go RIC I	Charger Case Ease	Charger Case Go	Combi RIC Charger V2	Combi BTE Charger V2
								
	075-0108-11090	075-0101-11	075-0057-11	075-0055-11	075-0080-11	075-3017-11	075-0063-11	075-0065-11
Phonak Audéo™ E-RSphere	•							
Phonak Audéo E-R	•							
Phonak CROS E-R	•							
Phonak Virto™ I-R		•						
Phonak Audéo I-RSphere			•					
Phonak Audéo I-R				•				
Phonak CROS I-R				•				
Phonak Audéo L-R/RT					•			
Phonak Audéo L-R Life						•		
Phonak Audéo L-R Fit					•			
Phonak CROS L-R					•			
Phonak Naída™ L-PR								•
Phonak Audéo P-R/RT							•	
Phonak Audéo P-R Life						•		
Phonak Audéo P-R Fit							•	
Phonak CROS P-R							•	
Phonak Naída P-PR								•
Phonak Audéo M-R/RT							•	
Phonak Bolero M-PR								•

● 071-0003 Power Pack compatible

Power Pack is not recommended to be charged with third-party battery chargers (power banks) and should only be charged with the combi case power supply. Power Pack always ships ground.

Charger compatibility guide.

	Charger Go RIC E	Charger Go Virto I	Charger Go RICSphere I	Charger Go RIC I	Charger Case Ease	Charger Case Go	Combi RIC Charger V2	Combi BTE Charger V2
								
	075-0108-11090	075-0101-11	075-0057-11	075-0055-11	075-0080-11	075-3017-11	075-0063-11	075-0065-11
External Power Supply (EPS) 5V US - no cable included 075-0100-11	•	•	•	•	•	•	•	•
USB A to USB C cable 043-3116	•	•	•	•	•			
USB microUSB 1m 28AWG No data 043-0002							•	•
USB A to USB C cable with ferrite bead 043-3116-1						•		

Power Supply

External Power Supply (EPS)
5V 1A US
075-0100-11



Cable USB Micro USB 1m
28AWG No data
043-0002



USB A to USB C with ferrite
043-3116-1



USB A to USB C
043-3116



028-6822-03/2026-10/hc © 2024-2026 Sonova AG. All rights reserved.